

## IL-1RL2 (C-Fc), Human, Recombinant

货号 : PCK147

### 产品信息

别名	Interleukin-1 Receptor-like 2; IL1RL2; IL-36 Receptor; IL-36R; IL-1Rrp2; IL1R-rp2;
物种	Human
表达宿主	Human Cells
序列信息	Asp20-Tyr337
检索号	Q9HB29
分子量	63.1 kDa
标签	C-Fc

### 产品特性

纯度	>95% as determined by reducing SDS-PAGE.
内毒素	<1.0 EU per $\mu\text{g}$ as determined by LAL test.
保存	Lyophilized protein should be stored at $-5\sim-20^{\circ}\text{C}$ , stable for one year after receipt. Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $-5\sim-20^{\circ}\text{C}$ for 3 months.
运输	Ambient temperature or ice pack.
制剂	Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of PBS, pH7.4.



## 复融

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## 背景介绍

Interleukin-1 Receptor-like 2 is a Protein that in humans is encoded by the IL1RL2 gene, belongs to the Interleukin-1 Receptor family. IL1RL2 is the Receptor for Interleukin-36 (IL36A, IL36B and IL36G). After binding to Interleukin-36 associates with the coReceptor IL1RAP to form the Interleukin-36 Receptor complex which mediates Interleukin-36-dependent activation of NF-kappa-B, MAPK and other pathways.

## SDS-PAGE

