

NT3/ NGF-2, Human, Recombinant

货号 : PCK210

产品信息

别名	Neurotrophin-3; NT-3; HDNF; Nerve Growth Factor 2; NGF-2; Neurotrophic Factor; NTF3
物种	Human
表达宿主	E.coli
序列信息	Tyr139-Thr257
检索号	P20783
分子量	13.6 kDa

产品特性

纯度	>95% as determined by reducing SDS-PAGE.
内毒素	<1.0 EU per µg as determined by LAL test.
保存	Lyophilized protein should be stored at -5~-20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -5~-20°C for 3 months.
运输	Ambient temperature or ice pack.
制剂	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 250mM NaCl, pH 7.2.



复融

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

背景介绍

Neurotrophin-3 (NT-3) is a member of the NGF family of neurotrophic factors and is structurally related to β -NGF, BDNF and NT-4. The NT3 cDNA encodes a 257 amino acid residue precursor Protein with a signal peptide and a pro Protein that are cleaved to yield the 119 amino acid residue mature NT3. The amino acid sequences of mature human, murine and rat NT-3 are identical. NT-3 selectively promotes the differentiation and survival of specific neuronal subpopulations in both the central as well as the peripheral nervous systems.

SDS-PAGE

